OPTIC & SCIENTIFIC TOYS

www.navir.it navir@navir.it

Via della Resistenza 34/B - 20090 Buccinasco (Milano) - Italy
Tel. 02-45700262 - Fax 02-4880736

KIT DISCOVERY

OPTICAL NATURE KIT

OPTIC & SCIENTIFIC TOYS

**Scientific consultancy and instruction booklet:
Cesare Baj – Cernobbio (Como) – Italy**

Graphic design: Mirko Alpiani (info@m-design.it)

ITEM 8093

MADE IN ITALY

ÓPTIC & SCIENTIFIC TOYS

Brief historical outline

Lenses were already known to the Greeks and Romans who used them to light fires (burning lenses). It was in 1300 that lenses were first produced in Italy to manufacture eyeglasses, a production that developed mainly in Holland.

It was precisely in Holland that in 1600 it was discovered that, by adequately combining together different lenses, it was possible to magnify either very distant or very small objects. This led to the creation of the microscope and the telescope, two instruments that have enabled humankind to in better know the world we live.

The microscope, first manufactured with a single lens (simple microscope) and later with two (compound microscope), has enabled us to study the hidden world which is not visible to the naked eye and, in particular, made finally visible formerly unknown forms of life.


Among the first to use the microscope were the Italian scientist Galileo Galilei (1564-1642) who actually gave the instrument its name (from the Greek mikros = little and skopein = to see), the Dutch Antony van Leeuwenhoek (1632-1723) and the English Robert Hooke (1635-1703). The latter discovered the existence of cells by observing a piece of cork oak bark.

Although Zoomscope may seem to be just a toy, it should be remembered that many important discoveries were made with instruments with similar power and characteristics.

Zoomscope

The Zoomscope consists of an eye lens, an eyepiece and a mobile lens that allows to vary the degree of magnifying power. The magnifying power can be varied by turning the wheel situated on the side of the instrument.

Usually one starts with the lowest magnification to increase it gradually, until an optimal view of the object is achieved.

Like all microscopes, Zoomscope produces an upside-down image of the observed object.

Specimen lighting

The Zoomscope allows two kinds of observation.

If the material is transparent, it can be observed in transmitted light, i.e. with the light that passes through the material. This kind of viewing is achieved by aiming the instrument towards a light source such as the sky or a lamp (1).

Vice versa, if the material is opaque, one should use incident light. In this case the built-in lighting system of the instrument will be switched on, thus lighting whatever has been placed within the field of vision (2).

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


When the instrument lighting system is not in use, it can be stored in its proper housing. The battery should be inserted as shown in the figure. To switch the bright spotlight bulb on and off, turn the small handle on the side opposite to the light.

CAUTION!

This functions with a battery LR1 (size N), 1,5 V.

Only a battery of the same type or equivalent to that recommended should be used. The rundown battery must be removed and thrown in the special differentiated waste container.

Focusing

Focusing is achieved by slowly rotating the transparent part (or, alternatively, by rotating the instrument).

If the object you are viewing is flat, it can be brought into focus with just one operation. If the object is thick, only one section at a time can be focused. Therefore, to see other sections of the object that are not on the same level, the focus must be adjusted continuously.

Slides

Sometimes the specimens are better observed on a slide or between two slides, which are supplied with the Zoomscope.

To observe a slide, the stage must be assembled to the instrument. The specimen is then applied to the slide (either dry or moistened with a drop of water). Finally, the slide is covered with another slide and the two are fixed together with adhesive tape.

The material can be better observed if it has previously been stained, for instance with iodine tincture.

The micrometer

One of the two supplied slides presents two small scales to measure the size of the observed object. One has a notch every half millimeter and the other has one every 32nd of an inch.


To measure the observed object, just place it next to the scale, which is used as if it were a ruler.

Viewing

The Zoomscope allows the observation of any object, live or inanimate, placed within the field of vision and brought into focus.

Hereunder we are listing a number of things that can be interesting to observe:

Caution!


Never aim the zoomscope (or any other optical instrument) towards the sun, for it can cause eye damage.

Insects

Insects offer many opportunities for observation. The one at the right is a small fly.


Leaves and other vegetable parts

Vegetables should be sliced as thin as possible. To this purpose professionals use an instrument called a microtome. You can use a sharp knife or a razor blade (always with the assistance of an adult).

Crystals

Dissolve a small amount of kitchen salt in little boiling water and stir for a long time. In this way you obtain a saturated salt solution. Apply one or two drops of this solution to the slide and allow all the water to evaporate. What is left after all the water has evaporated are the typical cubic crystals of sodium chloride, i.e. of common salt. You can follow the same procedure with sugar or any other crystalline substance.

Hair

Observe the difference between straight, curly or frizzy hair.

Stamps, notes, documents

The Zoomscope can be used to study and verify the authenticity of every kind of printed matter and documents.

Fingerprints

Police all over the world use the microscope to study and compare fingerprints.

Fabrics

Fabric weaves can be observed in detail.

Printed matter

Color prints are an interesting object to observe with the Zoomscope. With it you can learn how an infinite range of colors can be obtained simply by different combinations of the four basic colors: cyan, magenta, yellow and black.

Triple Bug Viewer

The instrument allows you to observe objects and small animals. The object can be observed from above, from the sides, through the transparent surfaces or from the bottom, through the mirror. By using the additional lenses, it is possible to increase the enlargement. If you place small animals onto the viewer, treat them with care, ensuring that you do not harm them. Free them immediately after the observation. The tweezers supplied are useful for handling small objects.

Megalens

Lenses were already used by the ancient Greeks and Romans to light fires (lighting lenses). They began to be produced for glasses around 1300, in Italy. Later, the glasses industry was developed in Holland and then in other European countries.

The Magnifying glass allows you to observe an enlarged object. It includes a second, smaller lens which allows you increase the enlargement. The pedestal has a pair of tweezers for handling small objects. By extracting the small pedestal and placing the magnifying glass onto a flat surface, it is possible to work on objects using both hands.

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